


visicolor® ECO Calcium Calcio


Messbereich / Range / Domaine de mesure / Rango

1 Δ \triangleq 5 mg/l Ca²⁺


REF 931 012




1
5 ml Probe
Sample
Echantillon
Muestra




2 Δ Ca-1




3
Umschwenken
Shake gently
Secouer légèrement
Mezclar volteándolo



4 Δ Ca-2
bis Farbumschlag
till coloration
jusque inversion de couleur
hasta cambio de color



5
1 Δ \triangleq 5 mg/l Ca²⁺
rot → blau
red → blue
rouge → bleu
rojo → azul



1+4
Meerwasser
Sea water
Eau de mer
Agua de mar