




Chromate / Cromato




Cat. No.	Cuvette	Measurement range	Method	Wavelength
918 25	10 mm	0.1 - 3.0 mg/l Cr (VI) 0.1 - 6.0 mg/l CrO ₄ ²⁻	(1) 251	540 nm
	50 mm	0.01 - 0.45 mg/l Cr (VI) 0.01 - 1.00 mg/l CrO ₄ ²⁻	(1) 251	




1 R1




2.0 ml R2




1'00 min
shake well
intense agiter
agitar intensamente




20 ml sample
échantillon
muestra




mix
mélanger
mezclar




fill up with dist. water
emplir avec l'eau dist.
llenar con agua dest.




mix
mélanger
mezclar




5'00 min




rinse and fill
rincer et emplir
lavar y llenar



measurement
mesure
medición



blank value
zéro
cero



blank measurement
zéro mesure
cero medición