




Hydrazine




Cat. No.	Cuvette	Measurement range	Method	Wavelength
918 44	10 mm	0.01 - 1.50 mg/l N ₂ H ₄	(1) 441	436 nm
	50 mm	0.002 - 0.250 mg/l N ₂ H ₄		




1
20 ml sample
échantillon
muestra




2
2.0 ml R1




3
mix
mélanger
mezclar




4
2.0 ml R2




5
mix
mélanger
mezclar




6
fill up with dist. water
emplir avec l'eau dist.
llenar con agua dest.




7
rinse and fill
rincer et emplir
lavar y llenar



6
10'00 min




7
measurement
mesure
medición



1 20 ml
dist. water
l'eau dist.
agua dest.
2 ... **6**
see above
cf ci-dessus
véase arriba

blank value
zéro
cero



7
blank measurement
zéro mesure
cero medición