

visicolor® ECO Ammonium 15

Amonio


Test 5-10

REF 931 210


585 nm

Method(e) / Método


5101	0.4–6.2 mg/L NH ₄ -N
5102	0.5–8.0 mg/L NH ₄ ⁺
5103	0.5–8.0 mg/L NH ₃




Probe / sample
échantillon / muestra




1
Küvetten spülen
Rinse test tubes
Rincer les éprouvettes
Lavar los tubos




2
5 mL Probe
Sample
Echantillon
Muestra



3
10 NH₄-1




4
Schütteln
Shake
Agiter
Agitar




5
1 NH₄-2



6
Schütteln
Shake
Agiter
Agitar




7
5'00 min




8
4 NH₄-3



9
Schütteln
Shake
Agiter
Agitar




10
Säubern
Clean
Nettoyer
Limpiar



11
7'00 min



12
Messung
Measurement
Mesure
Medición



1+9
Meerwasser
Sea water
Eau de mer
Agua de mar