

visocolor® ECO Chlor 2 (frei)

Chlorine (free) / Chlore (libre) / Cloro (libre)

Test 5-15/5-16


REF 931 215/931 216

540 nm


Method(e) / Método

5151 / 5161


0.10–2.00 mg/L Cl₂




Probe / sample
échantillon / muestra




1 3 Cl_2 -1




2 3 Cl_2 -2



3 5 mL Probe
Sample
Echantillon
Muestra



4 Schütteln
Shake
Agiter
Agitar



5 Säubern
Clean
Nettoyer
Limpiar



6 Messung sofort
Measurement immediately
Mesure immédiatement
Medición inmediatamente



Meerwasser
Sea water
Eau de mer
Agua de mar

visocolor[®] ECO Chlor 2 (gesamt)

Chlorine (total) / Chlore (totale) / Cloro (total)

Test 5-15


REF 931 215

540 nm


Method(e) / Método

5152


0.10–2.00 mg/L Cl₂




Probe / sample
échantillon / muestra




3 Cl₂-1




3 Cl₂-2



5 mL Probe
Sample
Echantillon
Muestra



Schütteln
Shake
Agiter
Agitar




3 Cl₂-3



Schütteln
Shake
Agiter
Agitar



Säubern
Clean
Nettoyer
Limpiar



2'00 min



Messung
Measurement
Mesure
Medición



Meerwasser
Sea water
Eau de mer
Agua de mar