

visocolor[®] ECO Chrom(VI)

Chromium(VI) / Chrome(VI) / Cromo(VI)

Test 5-20

REF 931 220

540 nm


Method(e) / Método

5201

0.02–0.50 mg/L Cr(VI)

5202


0.04–1.00 mg/L CrO₄²⁻




Probe / sample
échantillon / muestra



1
Küvetten spülen
Rinse test tubes
Rincer les éprouvettes
Lavar los tubos




2
5 mL Probe
Sample
Echantillon
Muestra



3
5  Cr-1




4
Schütteln
Shake
Agiter
Agitar




5
5  Cr-2



6
Schütteln
Shake
Agiter
Agitar




7
Säubern
Clean
Nettoyer
Limpiar



8
3'00 min



Messung
Measurement
Mesure
Medición



Meerwasser
Sea water
Eau de mer
Agua de mar