

visocolor[®] ECO Kupfer

Copper / Cuivre / Cobre

Test 5-37


REF 931 237

585 nm


Method(e) / Método

5371


0.1–5.0 mg/L Cu²⁺




Probe / sample
échantillon / muestra




1
Küvetten spülen
Rinse test tubes
Rincer les éprouvettes
Lavar los tubos




2
5 mL Probe
Sample
Echantillon
Muestra



3
5  Cu-1




4
Schütteln
Shake
Agiter
Agitar




5  Cu-2



6
Schütteln
Shake
Agiter
Agitar




7
Säubern
Clean
Nettoyer
Limpiar



8
10'00 min



9
Messung
Measurement
Mesure
Medición



Meerwasser
Sea water
Eau de mer
Agua de mar