

# visocolor<sup>®</sup> ECO Kalium

## Potassium / Potasio

# Test 5-32


REF 931 232

690 nm


Method(e) / Método

5321


2–25 mg/L K<sup>+</sup>




Probe / sample  
échantillon / muestra



**1**  
Küvetten spülen  
Rinse test tubes  
Rincer les éprouvettes  
Lavar los tubos




**2**  
10 mL Probe  
Sample  
Echantillon  
Muestra



**3**  
15  K-1




**4**  
Schütteln  
Shake  
Agiter  
Agitar



**5**  
1  K-2




**6**  
30 s Schütteln  
Shake  
Agiter  
Agitar



**7**  
Säubern  
Clean  
Nettoyer  
Limpiar



**8**  
Messung sofort  
Measurement immediately  
Mesure immédiatement  
Medición inmediatamente



1+1  
Meerwasser  
Sea water  
Eau de mer  
Agua de mar